



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® hinged lid jars 142 x 72 x 25 mm, 10  
pcs./pack**

**€27.20  
plus VAT &  
Shipping**

## Product Images



## Description

---

Practical storage boxes made of transparent, crystal clear polystyrene with hinged lid.  
Temperature resistant from -40°C to +70°C.

## Additional Information

No.	2-1512
Manufacturer (Brand)	neoLab
EAN	4058072057176
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	255 ml
sterile	No
Wide	142 mm
Form	rectangular
Height	25 mm
Length	72 mm
Depth	72 mm
Type closure/ lid	Slip closure
suitable for	variable applications
Cans and tins properties	with closure/lid
TBST MAX	70 °C
TBST MIN	-40 °C
Closure / lid properties	closed (attached)
suitable for cooling in	Freezer/ freezer
Color closure/ lid	transparent (crystal clear)
Material closure/ lid	Polystyrene (PS)

