



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® hinged lid jars 145 x 90 x 40 mm, 10
pcs./pack**

**€34.70
plus VAT &
Shipping**

Product Images



Description

Practical storage boxes made of transparent, crystal clear polystyrene with hinged lid.
Temperature resistant from -40°C to +70°C.

Additional Information

No.	2-1513
Manufacturer (Brand)	neoLab
EAN	4058072057183
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	522 ml
sterile	No
Wide	145 mm
Form	rectangular
Height	40 mm
Length	90 mm
Depth	90 mm
Type closure/ lid	Slip closure
suitable for	variable applications
Cans and tins properties	with closure/lid
TBST MAX	70 °C
TBST MIN	-40 °C
Color closure/ lid	transparent (crystal clear)
Material closure/ lid	Polystyrene (PS)
Closure / lid properties	closed (attached)
suitable for cooling in	Freezer/ freezer

