



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Usbeck jar with lid made of 18/8, 400 ml

€31.00
plus VAT &
Shipping

Product Images



Description

Stainless steel jars with button lids for storing sensitive samples or paste-like material. They are easy to clean, autoclavable and chemically resistant.

Material: 18/8 stainless steel

Additional Information

| | |
|--------------------------|---|
| No. | UB-0358 |
| Manufacturer (Brand) | Usbeck |
| old neoLab article no. | 1-1181 |
| EAN | 4058072401597 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| diameter | 85 mm |
| Color | silver |
| Material | Stainless steel |
| Volume | 400 ml |
| sterile | No |
| Form | around |
| Height | 85 mm |
| Type closure/ lid | Overlay closure |
| suitable for | variable applications |
| Cans and tins properties | with closure/lid food grade electrically conductive |
| Closure / lid properties | with handle / knob |
| suitable for cooling in | Refrigerator freezer/ freezer |
| Color closure/ lid | silver |
| Material closure/ lid | Stainless steel |

