




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Usbeck tin with lid made of 18/8, 1200 ml

**€46.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Stainless steel jars with button lids for storing sensitive samples or paste-like material. They are easy to clean, autoclavable and chemically resistant.

Material: 18/8 stainless steel

## Additional Information

---

No.	UB-0360
Manufacturer (Brand)	Usbeck
old neoLab article no.	1-1183
EAN	4058072401610
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	120 mm
Color	silver
Material	Stainless steel
Volume	1200 ml
sterile	No
Form	around
Height	120 mm
Type closure/ lid	Overlay closure
suitable for	variable applications
Cans and tins properties	with closure/lid food grade electrically conductive
Color closure/ lid	silver
Material closure/ lid	Stainless steel
Closure / lid properties	with handle / knob
suitable for cooling in	Refrigerator freezer/ freezer

