



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Usbeck jar with lid made of 18/8, 3000 ml

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Stainless steel jars with button lids for storing sensitive samples or paste-like material. They are easy to clean, autoclavable and chemically resistant.

Material: 18/8 stainless steel

Additional Information

No.	UB-0363
Manufacturer (Brand)	Usbeck
old neoLab article no.	1-1186
EAN	4058072401641
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	200 mm
Color	silver
Material	Stainless steel
Volume	2500 ml
sterile	No
Form	around
Height	100 mm
Type closure/ lid	Overlay closure
suitable for	variable applications
Cans and tins properties	with closure/lid food grade electrically conductive
Closure / lid properties	with handle / knob
suitable for cooling in	Refrigerator freezer/ freezer
Color closure/ lid	silver
Material closure/ lid	Stainless steel

