




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Round jars 35 x 13 mm, 10 pcs./pack

€7.39
plus VAT &
Shipping

Product Images



Description

Round jars with tight-fitting lid. For sample storage.
Temperature resistant: -40°C to +75°C.
Material: Polystyrene

Additional Information

No.	2-1500
Manufacturer (Brand)	neoLab
EAN	4058072057053
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	35 mm
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	12 ml
sterile	No
Form	around
Height	13 mm
Type closure/ lid	Overlay closure
suitable for	variable applications
Cans and tins properties	with special chemical resistance with closure/lid
TBST MAX	70 °C
TBST MIN	-20 °C
Closure / lid properties	closed
suitable for cooling in	Freezer/ freezer
Color closure/ lid	transparent (crystal clear)
Material closure/ lid	Polystyrene (PS)

