




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**neolab® Round jars 69 x 40 mm, 10 pcs./pack**

**€17.90**  
**plus VAT &**  
**Shipping**

**Product Images**

---

## Description

---

Round jars with tight-fitting lid. For sample storage.  
Temperature resistant: -40°C to +75°C.  
Material: Polystyrene

## Additional Information

---

No.	2-1506
Manufacturer (Brand)	neoLab
EAN	4058072057114
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	69 mm
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	150 ml
sterile	No
Form	around
Height	40 mm
Type closure/ lid	Overlay closure
suitable for	variable applications
Cans and tins properties	with special chemical resistance with closure/lid
TBST MAX	70 °C
TBST MIN	-20 °C
Closure / lid properties	closed
suitable for cooling in	Freezer/ freezer
Color closure/ lid	transparent (crystal clear)
Material closure/ lid	Polystyrene (PS)

