



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® glass jars with swing stopper, 1550 ml, 221 mm high, 2 pcs./pack

Special Price

**€6.00** was

**€9.00**

**plus VAT &  
Shipping**

### Product Images



## Description

---

Glass jars with glass lid and swing stopper with red rubber seal.  
The lid is firmly connected to the swing stopper.  
The firmly locking closure prevents unintentional opening of the jar.  
The rubber ring ensures good tightness.

## Additional Information

|                          |   |
|--------------------------|---|
| No.                      | 2-2268  |
| Manufacturer (Brand)     | neoLab  |
| EAN                      | 4058072060244   |
| Transport temperature    | Room temperature  |
| autoclavable at 121°C    | No  |
| Color                    | transparent (clear)   |
| Material                 | Soda-lime glass (normal glass) 1  |
| Volume                   | 1550 ml   |
| sterile                  | No  |
| Wide                     | 100 mm  |
| Form                     | square  |
| Height                   | 221 mm  |
| Length                   | 100 mm  |
| Depth                    | 100 mm  |
| Type closure/ lid        | Swing stopper   |
| suitable for             | variable applications   |
| Cans and tins properties | with special chemical resistance with closure/lid suitable for foodstuffs tight closing |
| Closure / lid properties | closed (attached)   |
| suitable for cooling in  | Refrigerator  |
| Color closure/ lid       | transparent (clear)   |
| Material closure/ lid    | Soda-lime glass (normal glass)  |

