



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Vitlab plastic jar w. hinged lid, 2 ml, 14 x 32 mm

€0.90 As low as
€0.69
**plus VAT &
Shipping**

Product Images



Description

round box made of PE, with attached snap-on lid

Additional Information

No.	LV-0918
Manufacturer (Brand)	Vitlab
old neoLab article no.	E-4260
EAN	4058072295981
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	14 mm
Color	transparent (milky)
Material	High-density polyethylene (HDPE)
Volume	2 ml
sterile	No
Form	around
Height	32 mm
Type closure/ lid	Snap lock
Cans and tins properties	with closure/lid tight closing
Closure / lid properties	closed (attached)
suitable for cooling in	Refrigerator
Color closure/ lid	transparent (milky)
Material closure/ lid	High-density polyethylene (HDPE)

