



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® jars 91 x 15 mm, contents 13 ml, 10 pcs./pack

**€12.50
plus VAT &
Shipping**

Product Images



Description

Airtight, leak-proof cans with attached lid. Firmly latching, tightly closing lid with 3-fold sealing system. For specimen collection, storage, transport and shipping. Comply with the following USA regulations: I.M.S.; F.D.A.; U.S.D.A.; U.S.P. Dimensions refer to inside height and top, inside diameter, respectively.

Additional Information

No.	2-5250
Manufacturer (Brand)	neoLab
EAN	4058072004514
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	15 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	13 ml
sterile	No
Form	around
Height	91 mm
Type closure/ lid	Snap lock
suitable for	variable applications
Cans and tins properties	with special chemical resistance with closure/lid suitable for foodstuffs HACCP certified
TBST MAX	110 °C
TBST MIN	-10 °C
Color closure/ lid	transparent (milky)
Material closure/ lid	High Density Polypropylene (HDPP)
Closure / lid properties	closed (attached)
suitable for cooling in	Freezer/ freezer

