



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® rectangular boxes 50 x 30 x 25 mm, 10  
pcs./pack**

**€8.22  
plus VAT &  
Shipping**

## Product Images



## Description

---

Rectangular jars made of crystal-clear polystyrene with tight-fitting lid.

For specimen storage.

Temperature resistant: -40°C to +70°C

## Additional Information

No.	2-1521
Manufacturer (Brand)	neoLab
EAN	4058072057237
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polystyrene (PS)
Volume	37 ml
sterile	No
Wide	50 mm
Width/inside	47 mm
Form	rectangular
Height	25 mm
Height inside	23.0 mm
Length	30 mm
Depth	30 mm
Type closure/ lid	Hinged latch/lid
suitable for	variable applications
Cans and tins properties	with closure/lid
Length inside	27 mm
TBST MAX	70 °C
TBST MIN	-40 °C
Closure / lid properties	closed (attached)
Material closure/ lid	Polystyrene (PS)

