



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® rectangular boxes 57 x 47 x 28 mm, 10
pcs./pack**

**€17.90
plus VAT &
Shipping**

Product Images



Description

Rectangular jars made of crystal-clear polystyrene with tight-fitting lid.

For specimen storage.

Temperature resistant: -40°C to +70°C

Additional Information

No.	2-1522
Manufacturer (Brand)	neoLab
EAN	4058072057244
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polystyrene (PS)
Volume	75 ml
sterile	No
Wide	57 mm
Width/inside	54 mm
Form	rectangular
Height	28 mm
Height inside	26.0 mm
Length	47 mm
Depth	47 mm
Type closure/ lid	Hinged latch/lid
suitable for	variable applications
Cans and tins properties	with closure/lid
Length inside	44 mm
TBST MAX	70 °C
TBST MIN	-40 °C
Material closure/ lid	Polystyrene (PS)
Closure / lid properties	closed (attached)

