



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Rectangular boxes 70 x 15 x 16 mm, 10 pcs./pack

€13.40
plus VAT &
Shipping

Product Images



Description

Rectangular jars made of crystal-clear polystyrene with tight-fitting lid.
For specimen storage.
Temperature resistant: -40°C to +70°C

Additional Information

No.	2-1524
Manufacturer (Brand)	neoLab
EAN	4058072057268
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polystyrene (PS)
Volume	17 ml
sterile	No
Wide	70 mm
Width/inside	67 mm
Form	rectangular
Height	16 mm
Height inside	14.0 mm
Length	15 mm
Depth	15 mm
Type closure/ lid	Hinged latch/lid
suitable for	variable applications
Cans and tins properties	with closure/lid
Length inside	12 mm
TBST MAX	70 °C
TBST MIN	-40 °C
Material closure/ lid	Polystyrene (PS)
Closure / lid properties	closed (attached)

