




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Rectangular jars 325 x 215 x 44 mm, for pouring out gels

**€30.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Rectangular jars made of crystal-clear polystyrene with tight-fitting lid.  
For specimen storage.  
Temperature resistant: -40°C to +70°C

## Additional Information

No.	2-1526
Manufacturer (Brand)	neoLab
EAN	4058072057275
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polystyrene (PS)
Volume	306 ml
sterile	No
Wide	325 mm
Width/inside	322 mm
Form	rectangular
Height	44 mm
Height inside	42.0 mm
Length	215 mm
Depth	215 mm
Type closure/ lid	Snap lock
suitable for	variable applications
Cans and tins properties	with closure/lid
Length inside	212 mm
TBST MAX	70 °C
TBST MIN	-40 °C
Material closure/ lid	Polystyrene (PS)
Closure / lid properties	closed (attached)

