



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Beta waste container round with bottle 2 l

€294.00
plus VAT &
Shipping

Product Images



Description

Transparent acrylic waste container with bottle.

Diameter 140 mm for 2 ltr bottle.

The wall thickness of the acrylic is 9.5 mm. The corresponding bottle made of PE for liquid waste is included. The bottle is placed in the container in such a way that it and its cap can be easily accessed. The profile of the lid prevents slipping and provides increased shock resistance. The container also accepts waste bags of size 20 x 30 cm without the bottle.

Additional Information

No.	CD-1787
Manufacturer (Brand)	Thermo Scientific
old neoLab article no.	4-1803
EAN	4066292480155
Transport temperature	Room temperature
diameter	140 mm
Color	transparent (clear)
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Form	cylindrical
suitable for	Bottles. round
Thickness	9.5 mm
for bottles with volume	2000 ml

