



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Cannula Waste Box 2000 ml

€7.63
plus VAT &
Shipping

Product Images



Description

Secure, hermetically sealed waste boxes for used, contaminated needles.

The universal pull-off/twist-off devices are suitable for all common types of needles and insulin syringes.

Also suitable for Luer-Lock syringes and scalpels.

Insertion opening Ø 50 mm.

The patented tongue ring allows the box to be closed easily, securely and liquid-tight.

The filling level is clearly visible.

Suitable for environmentally friendly incineration or autoclaving.

All boxes meet the requirements of TRBA/BGR 250.

TÜV/GS-tested.

Additional Information

No.	1-1543
Manufacturer (Brand)	neoLab
EAN	4058072032975
autoclavable at 121°C	Yes
diameter	112 mm
Color	yellow
Material	Polypropylene (PP)
Volume	2000 ml
Height	202 mm
Type closure/ lid	Pressure cap
suitable for	Syringes
Buckets and barrels properties	with closure/lid
DM top	140 mm
DM below	107 mm

