



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene vacuum balloon bottles for heavy duty, polypropylene 10 l Case

€112.80 As low
as
€91.33
plus VAT &
Shipping

Product Images

Description

Heavy Duty Nalgene™ Polypropylene Vacuum Balloon Bottle with Handle These heavy duty vacuum balloon bottles are suitable for the most extreme operating conditions, can be used as vacuum traps and hold a full vacuum for up to 8 hours. Capacity (Anglo-American): 2.6 gal. Closure Material: polypropylene Closure Size: 83 B Description: heavy duty vacuum balloon bottles Outer Diameter (metric): 250 mm Capacity (metric): 10 l Type: balloon bottle with handle For use with (application): Ideal for the most extreme operating conditions Contains: Handles Quantity per carton: 6 Quantity per package: 1 Volume to rim (Anglo-American): 3 gal. Volume to rim (metric): 12 l Autoclavable: Yes Unit Size: Case of 6

Additional Information

| | |
|-----------------------|--------------------|
| No. | CD-0532 |
| Manufacturer (Brand) | Nalgene |
| VGKL number | 228722261 |
| EAN | 4058072176952 |
| Gross weight | 9.48 |
| Transport temperature | Room temperature |
| Material | Polypropylene (PP) |
| Volume | 10000 ml |
| Country of origin | USA |
| Packaging volume | 0.18000000 cc |
| Customs tariff number | 39233090 |

