



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene autoclavable polypropylene balloon bottle with tap 50 l Each 83B

€293.00
plus VAT &
Shipping

Product Images

Description

Nalgene™ autoclavable polypropylene balloon bottle with stopcock, 10 l Thermo Scientific Nalgene autoclavable clarification bottles with stopcock have excellent chemical resistance and are suitable for distilled water as well as nutrient media. Capacity (Anglo-American): 13 gal . Closure Material: Polypropylene, Thermoplastic Elastomer Seal Closure Size: 83B Description: Autoclavable Balloon Bottle with Tap Inner Diameter (metric) Neck: 64 mm Type: Balloon Bottle with Tap For use with (Application): Long-term storage, and dispensing solutions and media Graduations: Permanent graduation in 5-liter or 1-gallon increments Height (metric) with cap: 678 mm Includes: Integrated shoulder handles, molded-in stopcock Number per carton: 1 Number per package: 1 Autoclavable: Yes Color: White (opaque) Unit Size: Each

Additional Information

No.	CD-0571
Manufacturer (Brand)	Nalgene
VGKL number	228723193
EAN	4058072981075
Gross weight	5.28
Transport temperature	Room temperature
Material	Polypropylene (PP)
Volume	50 liters
Country of origin	USA
Packaging volume	0.10300000 cc
Customs tariff number	39233090

