




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene round LDPE balloon bottles with tap 10 l Case of 6 83B

€98.60 As low
as
€79.33
plus VAT &
Shipping

Product Images



Description

Nalgene™ Round LDPE Balloon Bottles with Tap Use the Thermo Scientific™ Nalgene™ Round LDPE Balloon Bottles with Tap to collect and dispense distilled water, acids and reagents. Capacity (Anglo-American): 2.5 gal. Closure material: polypropylene, thermoplastic elastomer seal Closure size: 83 B Description: round balloon bottles with stopcock Inner diameter (metric) Neck: 64 mm Type: balloon bottle with stopcock For use with (application): For collecting and dispensing distilled water, acids and reagents Graduations: Permanent graduation in 5-liter or 1-gallon increments Height (metric) with cap: 389 mm Includes: Integral shoulder handles, molded-in stopcock Quantity per box: 6 Autoclavable: Not autoclavable Color: White (opaque) Unit Size: Case of 6

Additional Information

No.	CD-0565
Manufacturer (Brand)	Nalgene
VGKL number	228723182
EAN	4058072174521
Gross weight	12.4
Transport temperature	Room temperature
Material	Low density polyethylene (LDPE)
Volume	10000 ml
Country of origin	USA
Packaging volume	0.18200000 cc
Customs tariff number	39233090

