



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene round LDPE balloon bottles with tap 25 l Each 83B

€140.00
plus VAT &
Shipping

Product Images



Description

Nalgene™ Round LDPE Balloon Bottles with Tap Use the Thermo Scientific™ Nalgene™ Round LDPE Balloon Bottles with Tap to collect and dispense distilled water, acids and reagents. Capacity (Anglo-American): 6.5 gal. Closure material: polypropylene, thermoplastic elastomer seal Closure size: 83 B Description: round balloon bottles with stopcock Inner diameter (metric) Neck: 64 mm Type: Carboy with Spigot For use with (application): Collect and dispense distilled water, acids and reagents Graduations: Permanent graduation in 5-liter or 1-gallon increments Height (metric) with closure: 594 mm Includes: Integral shoulder handles, Molded-in spigot boss Autoclavable: Not autoclavable Color: White (opaque) Unit Size: Each

Additional Information

No.	CD-0567
Manufacturer (Brand)	Nalgene
VGKL number	228723184
EAN	4066292469365
Gross weight	2.9
Transport temperature	Room temperature
Material	Low density polyethylene (LDPE)
Volume	25 liters
Country of origin	USA
Packaging volume	0.05600000 cc
Customs tariff number	39233090

