



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene round LDPE balloon bottles with tap 4 l Case of 6 38 to 430 mm

€298.02
plus VAT &
Shipping

Product Images



Description

Nalgene™ Round LDPE Balloon Bottles with Tap Use the Thermo Scientific™ Nalgene™ Round LDPE Balloon Bottles with Tap to collect and dispense distilled water, acids and reagents. Capacity (Anglo-American): 1 gal. Cap Material: polypropylene, thermoplastic elastomer seal Cap Size: 38 to 430 mm Description: round balloon bottles with stopcock Inner Diameter (metric) Neck: 26 mm Outer Diameter (metric): 154 mm Capacity (metric): 4 l Type: Carboy with Spigot For use with (application): Collect and dispense distilled water, acids and reagents Graduations: Permanent graduation in 5-liter or 1-gallon increments Height (metric) with closure: 335 mm Includes: Integral shoulder handles, Molded-in spigot boss Number per box: 6 Autoclavable: Not Autoclavable Unit Size: Case of 6

Additional Information

No.	CD-0564
Manufacturer (Brand)	Nalgene
VGKL number	228723181
EAN	4066292469358
Gross weight	3.74
Transport temperature	Room temperature
Material	Low density polyethylene (LDPE)
Volume	4000 ml
Country of origin	Mexico
Packaging volume	0.05400000 cc
Customs tariff number	39233090

