



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene Autoclavable Balloon Bottles

€170.10 As low
as
€137.50
plus VAT &
Shipping

Product Images

Description

Nalgene™ Polypropylene, Balloon Bottle The Thermo Scientific™ Nalgene™ autoclavable balloon bottles are an excellent choice for large volumes of culture media, vaccine production, and storage of large volumes of media and chemicals. Capacity (Anglo-American): 5.3 gal. Cap Material: polypropylene Cap Size: 83 B Description: polypropylene balloon bottle with handle Outer Diameter (metric): 282 mm Capacity (metric): 20 l Type: balloon bottle with handle For use with (application): For use with large volumes of culture media, for vaccine production, for large volumes of media, and for chemical storage Height (metric) with cap: 526 mm Includes: Handles Quantity per package: 1 Volume to rim (anglo): 5.5 gal. Volume to rim (metric): 22.5 l Autoclavable: Yes Unit Size: Case of 4

Additional Information

No.	CD-0550
Manufacturer (Brand)	Nalgene
VGKL number	228722502
EAN	4066292469259
Gross weight	11.41
Transport temperature	Room temperature
Material	Polypropylene (PP)
Volume	20 liters
Country of origin	USA
Packaging volume	0.23800000 cc
Customs tariff number	39233090

