



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene Autoclavable Balloon Bottles

€100.40 As low
as
€81.17
plus VAT &
Shipping

Product Images

Description

Nalgene™ Polypropylene, Balloon Bottle The Thermo Scientific™ Nalgene™ autoclavable balloon bottles are an excellent choice for large volumes of nutrient media, vaccine production, and storage of large volumes of media and chemicals. Capacity (Anglo-American): 2.6 gal. Closure Material: polypropylene Closure Size: 83 B Description: polypropylene balloon bottle with handle Outer Diameter (metric): 250 mm Capacity (metric): 10 l Type: Carboy with Handle For use with (application): For use with large batches of culture media, vaccine product production, bulk media and chemical storage Height (metric) with closure: 389 mm Includes: Handles Quantity per package: 1 Volume to rim (Anglo-American): 3 gal. Volume to rim (metric): 12 l Autoclavable: Yes Unit Size: Case of 6

Additional Information

| | |
|-----------------------|--------------------|
| No. | CD-0549 |
| Manufacturer (Brand) | Nalgene |
| VGKL number | 228722501 |
| EAN | 4058072177645 |
| Gross weight | 11.25 |
| Transport temperature | Room temperature |
| Material | Polypropylene (PP) |
| Volume | 10000 ml |
| Country of origin | USA |
| Packaging volume | 0.18300000 cc |
| Customs tariff number | 39233090 |

