




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene Autoclavable Balloon Bottles

€249.00
plus VAT &
Shipping

Product Images

Description

Nalgene™ Polypropylene, Balloon Bottle The Thermo Scientific™ Nalgene™ autoclavable balloon bottles are an excellent choice for large volumes of nutrient media, vaccine production, and storage of large volumes of media and chemicals. Capacity (Anglo-American): 13 gal. Cap Material: polypropylene Cap Size: 83 B Description: polypropylene balloon bottle with handle Outer Diameter (metric): 379 mm Capacity (metric): 52 L Type: balloon bottle with handle For use with (application): For use with large volumes of culture media, for vaccine production, for large volumes of media, and for chemical storage Height (metric) with closure: 678 mm Includes: Handles Quantity per package: 1 Volume to rim (Anglo-American): 14 gal. Volume to rim (metric): 54 l Autoclavable: Yes Unit Size: Each

Additional Information

No.	CD-0551
Manufacturer (Brand)	Nalgene
VGKL number	228722503
EAN	4066292469266
Gross weight	4.63
Transport temperature	Room temperature
Material	Polypropylene (PP)
Volume	50 liters
Country of origin	USA
Packaging volume	0.06300000 ccm
Customs tariff number	39233090

