



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bürkle storage bottle without threaded neck, HDPE, 10 l, with cap

€46.35 As low
as
€39.70
plus VAT &
Shipping

Product Images

Description

Extra thick-walled, sturdy storage bottle with sturdy carrying handle for safe transport and storage. Printed with blue liter scale. The aeration/deaeration valve (accessory) allows air to enter the container to avoid a vacuum in the bottle. The canister can also be used very well as a clarification canister. A drain cock or other accessories with 3/4" female thread/ union nut can be connected via a threaded connection. HDPE bottle opening diameter 50 mm

Additional Information

No.	KL-2010
Dual use indicator	false
Manufacturer (Brand)	Bürkle
VGKL number	220530202
EAN	4058072278779
Transport temperature	Room temperature
Net weight	0.45650 kg
diameter	206 mm
Material	High density polyethylene (HDPE)
Volume	10000 ml
Country of origin	Germany
sterile	No
Customs tariff number	39233090
Height	427 mm

