



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Bucket made of PP, 12 l, blue

€27.80
plus VAT &
Shipping

Product Images



Description

For heavy duty industrial use!

FDA approved, therefore also suitable for the food industry.

The PP material remains hard and stable even when boiling water is used.

Temperature resistant from -20°C to +100°C.

With spout, unique pouring handle and stainless steel handle.

Measuring scale up to 12 l, total capacity 13 l.

Accessories: Precisely fitting lid made of PP, to protect the contents of the bucket (Order No. 2-7555 ff.).

Additional Information

No.	2-7553
Manufacturer (Brand)	neoLab
EAN	4058072009465
Transport temperature	Room temperature
diameter	325 mm
Color	blue
Material	Polypropylene (PP)
Volume	12 l
Form	around
Height	330 mm
Buckets and barrels properties	food safe cap/lid available as accessory
TBST MAX	100 °C
TBST MIN	-20 °C

