



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**neolab® Sterilizing pail with lid, stainless steel,  
outer dimensions 240 x 190 mm (ø x H)**

**€668.10  
plus VAT &  
Shipping**

## Product Images

---

## Description

---

The high-quality sterilization buckets with lids made of stainless steel are robust, practical and ideally suited for all autoclaving processes in the laboratory.

- available in different versions
- with inserted bottom
- 18/10 stainless steel
- Silk gloss brushed
- Perforation to be closed with lid
- Bucket handle made of flat steel
- Liquid tight

## Additional Information

---

No.	9-1200
Manufacturer (Brand)	neoLab
EAN	4058072991258
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	240 mm
Color	silver
Material	Stainless steel
Form	cylindrical
Height	190 mm
Type closure/ lid	Plug closure (inside)
TBST MAX	300 °C
Color closure/ lid	silver
Material closure/ lid	Stainless steel

