



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**neolab® sterilization pail with lid, stainless steel, outer dimensions 240 x 240 mm (ø x H)**

**€711.50  
plus VAT &  
Shipping**

## Product Images

---

## Description

---

The high-quality sterilization buckets with lids made of stainless steel are robust, practical and ideally suited for all autoclaving processes in the laboratory.

- available in different versions
- with inserted bottom
- 18/10 stainless steel
- Silk gloss brushed
- Perforation to be closed with lid
- Bucket handle made of flat steel
- Liquid tight

## Additional Information

---

|                       |                       |
|-----------------------|-----------------------|
| No.                   | 9-1201                |
| Manufacturer (Brand)  | neoLab                |
| EAN                   | 4058072991265         |
| Transport temperature | Room temperature      |
| autoclavable at 121°C | Yes                   |
| diameter              | 240 mm                |
| Color                 | silver                |
| Material              | Stainless steel       |
| Form                  | cylindrical           |
| Height                | 240 mm                |
| Type closure/ lid     | Plug closure (inside) |
| TBST MAX              | 300 °C                |
| Color closure/ lid    | silver                |
| Material closure/ lid | Stainless steel       |

