



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Sterilizing pail with lid, stainless steel,
outer dimensions 340 x 240 mm (ø x H)**

**€545.00
plus VAT &
Shipping**

Product Images

Description

The high-quality sterilization buckets with lids made of stainless steel are robust, practical and ideally suited for all autoclaving processes in the laboratory.

- available in different versions
- with inserted bottom
- 18/10 stainless steel
- Silk gloss brushed
- Perforation to be closed with lid
- Bucket handle made of flat steel
- Liquid tight

Additional Information

No.	9-1203
Manufacturer (Brand)	neoLab
EAN	4058072991289
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	340 mm
Color	silver
Material	Stainless steel
Form	cylindrical
Height	240 mm
Type closure/ lid	Plug closure (inside)
TBST MAX	300 °C
Color closure/ lid	silver
Material closure/ lid	Stainless steel

