




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bürkle filling aid for space saving canister 10l

€242.90
plus VAT &
Shipping

Product Images



Description

The filling aid makes it easy to fill liquids from Bürkle's 10 liter space-saving canister and dose them well. The space-saving canister sits on the indentations in the holder and is fixed with the tensioner. This makes it impossible for the canister to slip out. With the swiveling holder, the canister can be guided safely and the liquid medium can be filled and transferred without effort. The canister can be stored in the filling stand or replaced quickly and easily with the tensioner. The filling aid can be used wherever space-saving canisters are used, e.g. in laboratories and industry or in the food sector. Stainless steel V2A (1.4301) For Bürkle room-saving canisters 10 liters (PP, transparent or HDPE, black, antistatic) Scope of delivery: Filling aid (consisting of holder and filling stand), mountable plastic feet

Additional Information

| | |
|-----------------------|------------------|
| No. | KL-2226 |
| Manufacturer (Brand) | Bürkle |
| Dual use indicator | false |
| EAN | 4058072625092 |
| Transport temperature | Room temperature |
| Net weight | 2.69300 kg |
| sterile | No |
| Country of origin | Germany |
| Customs tariff number | 73269098 |

