



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN aluminum case unassembled for weight set stainless steel polished, 1 mg - 1 kg

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

---

- Aluminum case for weight sets with standard denomination of classes E1, E2
- Aluminum case for safe storage and transport in harsh industrial conditions

## Additional Information

---

No.	KP-1650
Manufacturer (Brand)	KERN
EAN	4058072128388
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
sterile	No
Weight max	2 kg
Weight min	1 kg
Manufacturer	CORE

