



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® wide-neck jug aluminum, 12.5 l, 218 mm neck Ø

€91.90
plus VAT &
Shipping

Product Images



Description

Aluminum transport cans made of 99.5% Al, according to EN 573.
Suitable for liquids, pastes, powders, granulates and much more.
Aluminum lid with stainless steel clamping ring.
Lid with inserted rubber seal.
Matt pickled outside and inside.

Additional Information

No.	7-3020
Manufacturer (Brand)	neoLab
EAN	4058072038502
Transport temperature	Room temperature
Color	silver
Material	Aluminum
Volume	12500 ml
sterile	No
Color bottle	silver
Height	375 mm
Type closure/ lid	Clamping ring closure
Jugs properties	with special thermal shock resistance with cap/lid with seal cap autoclavable bottle autoclavable
Neck size	far
Closure / lid properties	closed
DM bottle	220 mm
DM Neck	218 mm
Color closure/ lid	silver
Material bottle	Aluminum
Material closure/ lid	Aluminum

