



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® wide-neck jug aluminum, 30.5 l, 218 mm neck Ø

**€105.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Aluminum transport cans made of 99.5% Al, according to EN 573.  
Suitable for liquids, pastes, powders, granulates and much more.  
Aluminum lid with stainless steel clamping ring.  
Lid with inserted rubber seal.  
Matt pickled outside and inside.

## Additional Information

---

No.	7-3021
Manufacturer (Brand)	neoLab
EAN	4058072038519
Transport temperature	Room temperature
Color	silver
Material	Aluminum
Volume	30000 ml
sterile	No
Color bottle	silver
Height	480 mm
Type closure/ lid	Clamping ring closure
Jugs properties	with special thermal shock resistance with cap/lid with seal cap autoclavable bottle autoclavable
Neck size	far
Closure / lid properties	closed
DM bottle	300 mm
DM Neck	218 mm
Color closure/ lid	silver
Material bottle	Aluminum
Material closure/ lid	Aluminum

