



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Cooler Box IsoFreeze for 0°C, white

€165.40
plus VAT &
Shipping

Product Images



Description

For temperature-safe storage of PCR reagents, cells, enzymes or other temperature-sensitive materials.

For 20 reaction tubes 1.5/ 2.0 ml or

for 20 reaction tubes 0.5 ml (with inverted insert).

Cooling filling made of non-toxic gel.

Two models for different temperature ranges:

Order No. [2-9565](#)

white, for temperatures around 0°C,

keeps them for about 5 hours

Material: Polycarbonate

Additional Information

No.	2-9565
Manufacturer (Brand)	neoLab
EAN	4058072011178
Transport temperature	Room temperature
Color	White
Material	Polycarbonate (PC)
Wide	172 mm
Form	rectangular
Height	104 mm
Length	108 mm
Temperature holding time	0 °C for 5 hours
Depth	108 mm
Refrigerant medium	Refrigerant (closed system)
Cryogenic container features	with integrated refrigerant with closure/lid with grid insert insulating for different vessel types/volumes
Squares (WxL)	4 x 5
for vessels with DM up to	12.5 mm 5.5 mm
for vessels with volume up to	1.5 ml 2 ml 0.5 ml
for type (micro) vessel	PCR reaction vessels Reaction vessels
suitable for cooling in	Freezer/ freezer

