



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand Mini-Cooler 45 minutes for -70°C storage, 12x0.5 to 2.0 ml tubes

€109.65
plus VAT &
Shipping

Product Images



Description

For temperature control of a variety of reagents (enzymes, DNA, RNA, cell suspensions) by cold storage on the laboratory bench.

Made of durable polycarbonate, filled with non-toxic insulating gel.

Transparent lid.

Suitable for different storage temperature ranges.

For 12 tubes from 0.5 to 2.0 ml.

Order No. 2-2735

Color: white, 45 minutes for -70°C storage

Additional Information

No.	BN-0064
Manufacturer (Brand)	Brand
old neoLab article no.	2-2735
EAN	4058072196592
Transport temperature	Room temperature
Color	transparent (clear) white
Material	Polycarbonate (PC)
Wide	110 mm
Form	rectangular
Height	95 mm
Length	160 mm
Temperature holding time	-70 °C for 0.75 hours
Depth	160 mm
Refrigerant medium	Refrigerant (closed system)
Cryogenic container features	with integrated refrigerant with defined cooling rate with closure/lid with grid insert insulating for different vessel types/volumes alphanumerically coded
for vessels with volume up to	1 ml 1.2 ml 1.5 ml 2 ml 0.5 ml
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Squares (WxL)	3 x 4
for vessels with DM up to	13.0 mm 5.5 mm
for type (micro) vessel	Reaction vessels PCR reaction vessels

