



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® tabletop freezer box -10 to -50°C

**€4,469.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Small table-top units with low-noise refrigeration unit for the laboratory workplace.

- Interior free of internals, easy to clean
- Housing and interior made of stainless steel (4016, 4301)
- suitable for wall mounting
- optical and acoustic signal in case of deviation from the setpoint temperature
- Set point adjustment
- Rechargeable battery for approx. 60 hours failure indication
- Temperature recorder output 10 mV/k

### Technical data

Useful capacity: 35 liters

Outside dimensions:

580 x 765 x 540 mm (W x D x H)

Inside dimensions:

425 x 300 x 280 mm (W x D x H)

Weight: 81 kg

Ambient temp.: +12°C to +30°C

**Connection:** 230V/50Hz

**Order No. 7-4560**

**Temp. range:** -10°C to -50°C

## Additional Information

---

No.	7-4560
Manufacturer (Brand)	neoLab
EAN	4058072040628
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Temp MAX	-10 °C
Temp MIN	-50 °C
Volume	35 l
Wide	765 mm
Width/inside	300 mm
Height	540 mm
Height inside	280.0 mm
Length	580 mm
Depth	580 mm
Weight	81 kg
Length inside	425 mm
Voltage	220 V
for ambient temperature	12 to 30 °C
Energy supply	Power supply Battery
Functionality	electronic

