



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® tabletop freezer box -50 to -85°C

€5,776.00
plus VAT &
Shipping

Product Images



Description

Small table-top units with low-noise refrigeration unit for the laboratory workplace.

- Interior free of internals, easy to clean
- Housing and interior made of stainless steel (4016, 4301)
- suitable for wall mounting
- optical and acoustic signal in case of deviation from the setpoint temperature
- Set point adjustment
- Rechargeable battery for approx. 60 hours failure indication
- Temperature recorder output 10 mV/k

Technical data

Useful capacity: 35 liters

Outside dimensions:

580 x 765 x 540 mm (W x D x H)

Inside dimensions:

425 x 300 x 280 mm (W x D x H)

Weight: 81 kg

Ambient temp.: +12°C to +30°C

Connection: 230V/50Hz

Order no. 7-4562 **No. 7-4562**

Temp. range: -50°C to -85°C

Additional Information

No.	7-4562
Manufacturer (Brand)	neoLab
EAN	4058072040635
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Temp MAX	-50 °C
Temp MIN	-85 °C
Volume	35 l
Wide	765 mm
Width/inside	300 mm
Height	540 mm
Height inside	280.0 mm
Length	580 mm
Depth	580 mm
Weight	81 kg
Length inside	425 mm
Voltage	220 V
for ambient temperature	12 to 30 °C
Energy supply	Power supply Battery
Functionality	electronic

