



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Cooling Bag for Biological Samples, small

€48.40
plus VAT &
Shipping

Product Images



Description

These practical cooling bags are the ideal solution for transporting temperature-sensitive samples. The cooling function is activated with water. This activates a gel, which provides the cooling function. Filling up large Styrofoam boxes with ice is a thing of the past. Space and energy saving.

Simply submerge the inner bag completely in water for a maximum of 15 min and place the cooled samples inside. If the cooling wears off, simply moisten again. The outer pocket protects surrounding items from getting wet. Available in three different sizes and reusable.

Additional Information

No.	1-1404
Manufacturer (Brand)	neoLab
EAN	4058072090364
Color	purple
Wide	80 mm
Length	80 mm
Dimensions	80 x 180 mm
Type coolant	Cool bag

