



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**KGW Isotherm Dewar Vessel 9, approx. 1000 ml,
jacket blue, plastic-coated.**

**€168.18
plus VAT &
Shipping**

Product Images



Description

The Dewar vessels are made of borosilicate glass 3.3 (DIN/ISO 3585), silver-plated and evacuated. The Dewar vessels are supplied as standard in a blue, plastic-coated, galvanized sheet iron jacket. Ideal as cold trap and for all other applications where liquid nitrogen is required.

Type: 9; approx. 1000 ml

Additional Information

No.	GZ-0027
Manufacturer (Brand)	KGW Isotherm
old neoLab article no.	1-0503
EAN	4058072260941
Transport temperature	Room temperature
Color	blue
Material	Borosilicate glass steel
Volume	1000 ml
Height	290 mm
Height inside	235.0 mm
DM outside	105 mm
DM inside	77 mm
Cryogenic container features	insulating

