



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**KGW Isotherm Dewar Vessel 18, approx. 2500 ml,
jacket blue, plastic-coated.**

**€254.13
plus VAT &
Shipping**

Product Images



Description

The Dewar vessels are made of borosilicate glass 3.3 (DIN/ISO 3585), silver-plated and evacuated. The Dewar vessels are supplied as standard in a blue, plastic-coated, galvanized sheet iron jacket. Ideal as cold traps and for all other applications where liquid nitrogen is required.

Type 18; approx. 2500 ml

Additional Information

No.	GZ-0037
Manufacturer (Brand)	KGW Isotherm
old neolab article no.	1-0505
EAN	4058072187729
Transport temperature	Room temperature
Color	blue
Material	Borosilicate glass steel
Volume	2500 ml
Height	345 mm
Height inside	290.0 mm
DM outside	140 mm
DM inside	110 mm
Cryogenic container features	insulating

