



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**KGW Isotherm Dewar Vessel 3, approx. 500 ml,  
jacket blue, plastic-coated.**

**€127.65  
plus VAT &  
Shipping**

## Product Images



## Description

The Dewar vessels are made of borosilicate glass 3.3 (DIN/ISO 3585), silver-plated and evacuated. The Dewar vessels are supplied as standard in a blue, plastic-coated, galvanized sheet iron jacket. Ideal as cold trap and for all other applications where liquid nitrogen is required.

**Type: 3; approx. 500 ml**

## Additional Information

No.	GZ-0044
Manufacturer (Brand)	KGW Isotherm
old neoLab article no.	1-0501
EAN	4058072260965
Transport temperature	Room temperature
Color	blue
Material	Borosilicate glass steel
Volume	500 ml
Height	260 mm
Height inside	210.0 mm
DM outside	81 mm
DM inside	57 mm
Cryogenic container features	insulating

