



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**KGW Isotherm Dewar Vessel 6, approx. 800 ml,  
jacket blue, plastic-coated.**

**€142.00  
plus VAT &  
Shipping**

## Product Images



## Description

The Dewar vessels are made of borosilicate glass 3.3 (DIN/ISO 3585), silver-plated and evacuated. The Dewar vessels are supplied as standard in a blue, plastic-coated, galvanized sheet iron jacket. Ideal as cold traps and for all other applications where liquid nitrogen is required.

**Type: 6; approx. 800 ml**

## Additional Information

No.	GZ-0047
Manufacturer (Brand)	KGW Isotherm
old neoLab article no.	1-0502
EAN	4058072260972
Transport temperature	Room temperature
Color	blue
Material	Borosilicate glass steel
Volume	800 ml
Height	295 mm
Height inside	240.0 mm
DM outside	90 mm
DM inside	67 mm
Cryogenic container features	insulating

