



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene Dewar jar 1 l, inner dchm. 9.5 cm, inner height 19.4 cm

€290.40 As low
as
€235.25
plus VAT &
Shipping

Product Images



Description

Plastic Dewar flask, suitable for cryogenic temperatures.

Unbreakable, absolutely safe for handling liquid nitrogen, dry ice, etc.

Temperature resistant from -196°C to +100°C.

Double-walled filled with highly insulating PU foam, with ventilated insulated lid and carrying handle

Storage times for liquid nitrogen:

1 l = 7 h

Additional Information

No.	CD-1462
Manufacturer (Brand)	Nalgene
old neoLab article no.	1-0040
EAN	4058072219291
Transport temperature	Room temperature
Color	beige
Material	High-density polyethylene (HDPE)
Volume	1000 ml
Height	251 mm
Height inside	194.0 mm
Type closure/ lid	Overlay closure
Cryogenic container features	with closure/lid insulating
DM outside	200 mm
DM inside	95 mm
suitable for cooling in	Nitrogen tank

