



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Nalgene Dewar jar 4 l, inner dchm. 14.6 cm, inner height 28.7 cm

**€338.40** As low  
as  
**€274.00**  
plus VAT &  
Shipping

### Product Images



## Description

---

Plastic Dewar flask, suitable for cryogenic temperatures.

Unbreakable, absolutely safe for handling liquid nitrogen, dry ice, etc.

Temperature resistant from -196°C to +100°C.

Double-walled filled with highly insulating PU foam, with ventilated insulated lid and carrying handle

Storage times for liquid nitrogen:

4 l = 9 h

## Additional Information

---

No.	CD-1464
Manufacturer (Brand)	Nalgene
old neoLab article no.	1-0042
EAN	4058072219314
Transport temperature	Room temperature
Color	beige
Material	High-density polyethylene (HDPE)
Volume	4000 ml
Height	353 mm
Height inside	287.0 mm
Type closure/ lid	Overlay closure
DM outside	250 mm
DM inside	146 mm
Cryogenic container features	with closure/lid insulating
suitable for cooling in	Nitrogen tank

