



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Anesthesia/Induction Chamber for Mice

**€449.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Anesthesia box for mice:

High-quality plexiglass box glued water- and gas-tight. The supply of the anesthetic gas takes place in the lower area of the box, the suction takes place from the top. The hose connections are fixed as standard for hoses with an inner diameter of 6 mm. The hose connection size can be changed, as a 1/8" thread is used.

The cover is designed as a sliding cover and can be pulled out completely. Due to the sliding lid function, a minimum of anesthetic gas is swirled and thus only minimally brought to the outside. The sliding lid has a handle for easy sliding of the lid.

Dimension anesthesia box mouse inside (LxWxH): 153x78x100 mm

## Additional Information

---

No.	6-2451
Manufacturer (Brand)	neoLab
EAN	4058072161101
Transport temperature	Room temperature
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Height	105 mm
Height inside	100.0 mm

