



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Marienfeld slide storage box f. 12 slides, 100 x 85 x 35 mm

€14.40 As low
as
€11.00
plus VAT &
Shipping

Product Images



Description

Sturdy design with wooden frame and black cover, metal hinges and reliable snap closure.
Equally suitable for storage or transport.
Plastic holding rack and table of contents in lid for identification of slides 76 x 26 mm.

Additional Information

No.	MD-1530
Manufacturer (Brand)	Marienfeld
old neoLab article no.	1-6255
EAN	4058072320874
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	Black
Material	Wood brass
Wide	85 mm
Height	35 mm
Length	100 mm
Depth	100 mm
Type closure/ lid	Hinged latch/lid
suitable for	Microscope slides (76 x 26mm)
Boxes properties	stackable/collapsible with closure clip with labeling field suitable for shipping coded
Type Coding	indexable
for slide dimensions	76 x 26 mm

