



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Marienfeld Slide Storage Box f. 25 slides, 140 x 100 x 35 mm

**€19.80** As low as  
**€14.70**  
plus VAT &  
Shipping

### Product Images



## Description

---

Sturdy design with wooden frame and black cover, metal hinges and reliable snap closure.  
Equally suitable for storage or transport.  
Plastic holding rack and table of contents in lid for identification of slides 76 x 26 mm.

## Additional Information

---

No.	MD-1531
Manufacturer (Brand)	Marienfeld
old neoLab article no.	1-6257
EAN	4058072320881
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	Black
Material	Wood brass
Wide	100 mm
Height	35 mm
Length	140 mm
Depth	140 mm
Type closure/ lid	Hinged latch/lid
suitable for	Microscope slides (76 x 26mm)
Boxes properties	stackable/collapsible with closure clip with labeling field suitable for shipping coded
Type Coding	indexable
for slide dimensions	76 x 26 mm

