




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine Slide Box for 100 OT, ABS, yellow

**€10.90** As low as  
**€9.25**  
plus VAT &  
Shipping

### Product Images



## Description

---

Slide boxes for 100 slides 76 x 26 mm, with locking clip.  
No damage due to cork insert at the bottom, with hinged lid.  
Index directory in lid , stackable.

## Additional Information

---

No.	2-2437
Manufacturer (Brand)	neoLabLine
EAN	4058072060992
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	yellow
Material	Polystyrene (PS)
Wide	179 mm
Height	35 mm
Length	210 mm
Depth	210 mm
Type closure/ lid	Hinged latch/lid
suitable for	Microscope slides (76 x 26mm)
Boxes properties	stackable/collapsible with locking clip with labeling field coded
Type Coding	indexable
for slide dimensions	76 x 26 mm

