




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine Slide Box for 100 OT, ABS, blue

€10.90 As low as
€9.25
plus VAT &
Shipping

Product Images



Description

Slide boxes for 100 slides 76 x 26 mm, with locking clip.
No damage due to cork insert at the bottom, with hinged lid.
Index directory in lid , stackable.

Additional Information

No.	2-2438
Manufacturer (Brand)	neoLabLine
EAN	4058072061005
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	blue
Material	Polystyrene (PS)
Wide	179 mm
Height	35 mm
Length	210 mm
Depth	210 mm
Type closure/ lid	Hinged latch/lid
suitable for	Microscope slides (76 x 26mm)
Boxes properties	stackable/collapsible with locking clip with labeling field coded
Type Coding	indexable
for slide dimensions	76 x 26 mm

