



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLabLine Slide Box for 50 pcs., PS, red

€9.01 As low as
€8.20
plus VAT &
Shipping

Product Images



Description

Polystyrene slide boxes for 50 slides 76 x 26 mm.

With tight-fitting slip lid, for storage or shipping.

With numbering on the bottom for identification of the individual slides.

Stackable.

External dimensions: 177 x 83 x 32 mm (L x W x H)

Suitable for storage in LaboBox® magazine 2-2318.

Additional Information

No.	2-2348
Manufacturer (Brand)	neoLabLine
EAN	4058072060688
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	red
Material	Polystyrene (PS)
Wide	83 mm
Height	32 mm
Length	177 mm
Depth	177 mm
Type closure/ lid	Slip closure
suitable for	Microscope slides (76 x 26mm)
Boxes properties	stackable/collapsible suitable for shipping coded
Type Coding	numbered
for slide dimensions	76 x 26 mm

